

BIEBBING CONTROL KITS

BE PREPARED. PRESERVE LIFE.

Severe blood loss, often caused by motor vehicle accidents, acts of violence and natural disasters, can result in death within minutes. The SAM Medical Bleeding Control Kit includes essential medical devices designed to limit blood loss in emergency and traumatic situations. After calling emergency personnel, this life-preserving kit may help stabilize the injured in such events.



FEATURED TOURNIQUET

SAM XT EXTREMITY TOURNIQUET - DoD STOP THE BLEED APPROVED

Engineered for rapid application, SAM XT's innovative design normalizes the number of windlass turns needed compared to traditional tourniquets, enabling an easier, faster training of the product and intuitive use.





SAM[®]Bleeding Control Kit

KIT CONTENTS:

SAM XT EXTREMITY TOURNIQUET

Slack is the main cause of failed tourniquet application. SAM XT's TRUFORCE™ Buckle technology auto-locks at a predetermined amount of circumferential force, eliminating nearly all tourniquet slack. Simply click, twist, and secure for a life-preserving tourniquet application. The SAM XT is on the list of approved articles by the DoD for the Stop The Bleed Initiative.

MEDICAL GLOVES

Nitrile protective medical gloves - Latex free, powder free, size L.

MEDICAL SHEARS

Stainless steel medical shears allow the first responder to expose the wound and provide for appropriate treatment more efficiently.

COMPRESSED GAUZE

Designed to treat wounds involving moderate to severe blood loss. Can be secured in place with Emergency Bandage. 100% cotton, 6-ply, sterile.

5 EMERGENCY BANDAGE

Pressure bandage can be used directly over a wound or to hold gauze dressing in place with or without a pressure bar.

6 INSTRUCTION CARD

Includes detailed instructions for use and a diagram for proper tourniquet placement.







PRODUCT SPECIFICATIONS

PACKAGED SIZE & WEIGHT

Length: 9.7 in (24.5 cm) Height: 7 in (18 cm) Width: 3.2 in (8 cm) Weight: 10.9 oz (310 g)





SAM Medical Products® 12200 SW Tualatin Rd., Ste 200 Tualatin, OR 97062 USA

OPERATING SPECIFICATIONS

After calling emergency personnel, this life-preserving kit may help stabilize the injured in traumatic events.

Vacuum packed - KT901-V-EN U.S. Toll Free: 800.818.47 Non-Vacuum packed - KT901-NV-EN U.S. Local: 503.639.5474

GLOBAL SALES

U.S. Toll Free: 800.818.4726 U.S. Local: 503.639.5474 U.S. Fax: 503.639.5425 sammedical.com

 Size and weight vary by configuration and manufacturing process. ** Patent Pending SAM® is a registered trademark in the United States and other countries.
©2019 SAM Medical Products. All rights reserved.