



ChitoSAM™

[kahy-tuh-sam]

COVER. CLOT. CONTROL.

STOP BLEEDING RAPIDLY


ChitoSAM™ 100 is a high-performance hemostatic, non-woven chitosan dressing spun directly from chitosan derived crustaceans or snow crab shells. Designed to stop lethal bleeding rapidly, its ease of use is extremely effective. Works independently from the body's normal clotting processes.



FEATURED TECHNOLOGY

TEAR ANYWHERE PACKAGING

Innovative package design allows for quick access to ChitoSAM™ 100, tearable from any side.





ENGINEERED TO PRESERVE LIFE

- 1 PERFORMANCE GAUZE**
 ChitoSAM™ 100 has a hemostatic performance equivalent to products impregnated or coated with chitosan hemostatic agents.
- 2 CHITOSAN HEMOSTATIC AGENT**
 Unlike most hemostatic dressings that are only coated or impregnated with chitosan, ChitoSAM™ 100 is made from non-woven chitosan gauze, spun directly from chitosan derived from crustaceans or snow crab shells.
- 3 TEAR ANYWHERE PACKAGING**
 Innovative package design allows for quick access to ChitoSAM™ 100, tearable from any side.
- 4 Z-FOLD**
 Available for large wounds in 3in x 4ft and 3in x 6ft.
- 5 SINGLE PLY**
 Smaller 4in x 4in single ply ChitoSAM is able to be placed on wounds in addition to trauma gauzes.



PRODUCT SPECIFICATIONS

CHITOSAM™ 100 MODELS AVAILABLE (NON-CE)

- ChitoSAM 100 4 in x 4 in (10 cm x 10 cm) (CT100-A-EN)
- ChitoSAM 100 3 ft x 4 ft (7.6 cm x 1.22 m) (CT100-B-EN)
- ChitoSAM 100 3 ft x 6 ft (7.6 cm x 1.83 m) (CT100-C-EN)

Manufactured by SAM Medical Products
 12200 SW Tualatin Rd. Ste 200, Tualatin, OR, 97062 USA

CHITOSAM™ 100 MODELS AVAILABLE (CE)

- ChitoSAM 100 4 in x 4 in (10 cm x 10 cm) (CT101-A-EN)
- ChitoSAM 100 3 ft x 4 ft (7.6 cm x 1.22 m) (CT101-B-EN)
- ChitoSAM 100 3 ft x 6 ft (7.6 cm x 1.83 m) (CT101-C-EN)

Manufactured by BenQ Materials Corporation
 29, Jianguo E. Road, Gueishan
 33341, Taoyuan, Taiwan, R.O.C. +886.3.374.8800

GLOBAL SALES

U.S. Toll Free: 800.818.4726
 U.S. Local: 503.639.5474
 U.S. Fax: 503.639.5425
sammedical.com

WARNINGS

Do not place over eyes.
 Contains Chitosan from shellfish.

		STERILE R
		SAM® is a registered trademark in the United States and other countries. ©2019 SAM Medical Products. All rights reserved.



MedNet GmbH, Borkstrasse 10,
 48163 Muenster, Germany